

# IHE Integration Statement

Vendor	Product Name	Version	Date
Intuitive Imaging Informatics, LLC	ImageQube PACS	9.0	March 5, 2008
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Manager Image Creator Image Display	None
Patient Information Reconciliation	Image Manager	None
Consistent Presentation of Images	Image Manager Image Creator Image Display	None
Presentation of Grouped Procedures	Image Manager Image Display	None
Access to Radiology Information	Image Manager Image Display	None
Key Image Notes	Image Manager Image Creator Image Display	None
Portable Data for Imaging	Portable Media Creator Portable Media Importer Image Display	None

## Link to Standards Conformance Statements

<b>IHE</b>	<a href="http://www.IntuitiveImaging.com/support/compliance/I3IHEIntegrationStatement.pdf">www.IntuitiveImaging.com/support/compliance/I3IHEIntegrationStatement.pdf</a>
<b>HL7</b>	<a href="http://www.IntuitiveImaging.com/support/compliance/I3HL7ConformanceStatement.pdf">www.IntuitiveImaging.com/support/compliance/I3HL7ConformanceStatement.pdf</a>
<b>DICOM</b>	<a href="http://www.IntuitiveImaging.com/support/compliance/I3DicomConformanceStatement.pdf">www.IntuitiveImaging.com/support/compliance/I3DicomConformanceStatement.pdf</a>

## Link to general information on IHE

In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>
--	---	--

# Glossary

**Actors:**

Information systems or components of information systems that produce, manage, or act on categories of information required by operational activities in the enterprise.

**DICOM:**

Digital Imaging and Communications in Medicine. The industry standard for transfer of radiologic images and other medical information. DICOM enables digital communication between diagnostic and therapeutic equipment and systems from various manufacturers.

**HL-7:**

Health Level 7. A communications standard that enables healthcare applications to exchange key sets of clinical data. Health Level 7 is one of several ANSI accredited standards operating in the healthcare arena.

**IHE:**

Integrating the Healthcare Enterprise. A set of standards to assist healthcare enterprises in achieving systems integration, enabling secure access to vital information for optimal patient care. Sponsored jointly by the Radiological Society of North America (RSNA) and the Healthcare Information and Management Systems Society (HIMSS).

**Connectathon:**

The annual IHE event at which systems are tested.

**Transactions:**

Interactions between actors that communicate the required information through standards based messages.